Innovative Powder Strength Tester



Sample Size	Consolidation Pressure	Quick Results
The old way: 50 to 600 grams per test	The old way: 1000 Pa (1kPa) or more	The old way: as many as 5 tests to get one point
The new way: 0.06 to 0.1 gram per test	The new way: from 40 Pa to 50 kPa	The new way: one test = one point



Comes Complete With

- pre-programmed laptop
- bonus testing cell
- instruction manual
- one-year limited warranty on parts and labor
- micro-spatula for easy fill
- CE compliant



Manufactured and Distributed by:

Material Flow Equipment, LLC 5921 N County Road 225 Gainesville, FL 32609 USA Direct: 352-303-9123